BEC SARS-CoV-2 RT-LAMP kit for environmental samples

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BEC SARS-CoV-2 RT-LAMP kit
for environmental samples
P05302.10 ea

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The test uses licenses related to patents WO 00/28082, WO 01/34790, WO 01/77317, WO 02/24902 .

Fabriqué en France
Made in France

Référence : P11302

Version: 0120

Summary - Sommaire

▶ ENGLISH	4
• INTENDED USE.....	5
• CONTENT OF THE KIT	5
• MATERIALS REQUIRED BUT NOT PROVIDED	6
• WARNINGS AND PRECAUTIONS	7
• GENERAL INFORMATION	8
• PRODUCT DESCRIPTION	9
• OPERATING PROTOCOL	12
• INTERPRETATION OF RESULTS.....	24
• OPTIONAL: EXPORTING RESULTS	26
• TECHNICAL SUPPORT	28
• BIBLIOGRAPHY	29
▶ FRANÇAIS	31
• USAGE PREVU.....	32
• COMPOSANTS DU KIT	32
• MATERIEL REQUIS NON FOURNI.....	33
• MISES EN GARDE ET PRECAUTIONS.....	34
• INFORMATIONS GENERALES	36
• DESCRIPTION DU PRODUIT.....	37
• MODE OPÉRATOIRE	40
• INTERPRÉTATION DES RÉSULTATS	52
• OPTION : EXPORT DES RÉSULTATS	54
• ASSISTANCE TECHNIQUE	56
• BIBLIOGRAPHIE	57

English

► **Intended Use**

The BEC SARS-CoV-2 RT-LAMP kit for environmental samples is a test for rapid detection of SARS-CoV-2 combining a technology of rapid extraction of RNA from environmental samples and a RT-LAMP molecular test targeting 2 genes of SARS-CoV-2.

The analysis of this test must be carried out with the BEC Reader device.

BEC SARS-CoV-2 RT-LAMP kit for environmental samples has been developed for a use with air samples collected with Coriolis® air samplers and with surface samples.

► **Content of the kit**

Number of tests in a kit : 10.

Kit storage temperature : +4°C.

Kit shelf life : expiry date indicated on the outside of the box.

Kit content :

- 5 individual pouches containing 2 tests each.
- A cardboard rack, which makes it easier to handle the tubes, see Figure 2.
- A user guide.

Content of a pouch :

Description	Reference	Color	Quantity per pouch	Volume
BEC SARS-CoV-2 Assay Tubes for environmental samples	P04302	-	1 strip of 6 microtubes	-
BEC SARS-CoV-2 Lysing Buffer (green)	P13300	Green	2	1,5 mL
BEC SARS-CoV-2 Neutralization Buffer (blue)	P20300	Blue	2	1,5 mL
BEC SARS-CoV-2 Reaction Buffer (orange)	P21300	Orange	2	150 µL
BEC SARS-CoV-2 Transfer Tubes (violet)	P30300	Violet	2	-
sterile 3mL syringe	P31300	-	2	-
sterile needle for 3mL syringe	P32300	-	2	-
5 micron filter for 3mL syringe	P33300	-	2	-

► Materials required but not provided

- The BEC field reader pack including :
 - o The BEC Reader device (see Picture 2)
 - o A 25 µL fixed-volume pipette (orange)
 - o A 100 µL fixed-volume pipette (violet)
 - o A box of disposable filter-tips (recommended CLEARLINE cat. number 713116)
 - o Biological waste bin

- Sampling materials :
 - o Swabs for surface sampling
 - o Coriolis® air sampler and consumable (Bertin Technologies)
- Personal protective equipment for the user

► **Warnings and precautions**

Please read the user manual carefully before using the product.

- Before any use, please check that the product and its contents are not deteriorated and that all the contents are in the kit.
- Do not use a test from a pouch that has been opened for more than 24 hours, regardless of the storage temperature.
- Handle the samples as if they were infectious and/or dangerous, in accordance with the security procedures applicable in the laboratory.
- Wear adequate personal protective equipment (laboratory coat, disposable powder-free gloves, mask).
- Avoid microbiological and amplicon contaminations (by DNase/ RNase) of the sample and of the kit contents.
- Always use pipettes with disposable DNase/RNase-free filter-tips
- Always wear protective powder-free gloves when handling the contents of the kit.
- Do not open the reaction tubes after amplification to avoid any contamination with amplicons.
- Do not use the kit content after their expiry date.
- Dispose of the samples, contents in the pouch and waste material resulting from the test in

accordance with local safety regulations.

- This kit is intended for use with the reader approved for this purpose. Any directly interpretable result is dependent on the reader settings. If another reading instrument is used, it is the responsibility of the operator to ensure the correct calibration of that instrument. Bertin Technologies shall bear no responsibility for results obtained with any other reference system.
- Do not leave the microtubes in the reader once the reaction is completed.
- It is recommended to use a flat, clean and clear working surface in a covered and closed area at a temperature between 15°C and 25°C.
- Before opening the preparation tubes, make sure that all the liquid or lyophilized pellet contained in the tube is at the bottom of the tube. If needed, tap the tube several times downwards to lower any drops.

Never reopen the microtubes once they have been inserted in the BEC Reader device, nor after the end of the reaction. Re-opening under these conditions would result in environmental contamination. This contamination may cause false positive results for subsequent tests.

► General information

The BEC SARS-CoV-2 RT-LAMP kit for environmental samples is a field molecular test kit allowing the detection of the SARS-CoV-2 virus in environmental samples such as air or surface samples.

The SARS-CoV-2 is a coronavirus responsible for the emergence of a new respiratory infectious disease called Covid-19 that has caused a global pandemic. This disease, which appeared in China at the end of 2019, is transmitted mainly by air, with a high reproduction rate of between 2 and 3 in the absence of control and prevention measures. The severity of the clinical signs of the disease varies greatly from one individual to another, with a preponderance of severe forms in elderly individuals or those suffering from other diseases.

The scientific data available to date make it possible to estimate that 30% to 60% of infected subjects are asymptomatic but can transmit the disease. In this context, it is essential to be able to carry out diagnostic tests on a massive scale in order to limit the spread of the virus.

The BEC-SARS-CoV-2 kit has been developed in order to allow molecular tests to be carried out by non-specialized personnel directly on the field.

Source : Covid-19 Disease Sheets – Institut Pasteur – May 2020 (<https://www.pasteur.fr/en/medical-center/disease-sheets/covid-19-disease-novel-coronavirus>)

► **Product description**

The BEC SARS-CoV-2 RT-LAMP kit for environmental samples is an ensemble of reagents allowing rapid extraction of viral RNA and real-time RT-LAMP amplification of 2 sequences of the SARS-CoV-2 virus genome. It contains an exogenous positive control to ensure the correct running of the amplification.

RT-LAMP is a method of specific isothermal amplification of a nucleotide sequence with incorporation of a fluorescent nucleic acid stain during amplification (Li & Macdonald, 2015; Nagamine et al., 2002; Notomi et al., 2000; Yan et al., 2014). It amplifies a target sequence with two to three sets of primers at a constant temperature of 65°C. Compared to traditional detection methods, RT-LAMP is more specific and faster (Aebischer et al., 2014). It produces more amplicons and is less sensitive to inhibitors that are naturally present in biological matrices (Chander et al., 2014; Francois et al., 2011; Nie et al., 2012).

The BEC SARS-CoV-2 RT-LAMP kit for environmental samples allows 2 samples to be tested in parallel.

An analysis with the BEC SARS-CoV-2 RT-LAMP kit for environmental samples consists of carrying out the following steps :

- For the 2 samples in a sequential manner :
 - o Rapid extraction of viral RNA with chemical lysis and precipitation by chaotropic agents.
 - o Neutralization of the amplification inhibitors present in the lysate.
 - o Filtration of the lysate to remove aggregates of cell debris and reagents precipitated during neutralization.
- For the 2 samples in parallel: RT-LAMP amplification with the reaction buffer and the primers distributed in the microtubes (assay tubes), see Figure 1.


Figure 1: A strip of 6 microtubes allowing the amplification of 3 sequences of alien (C+ : positive control), RdRP and N genes for 2 samples

The test allows the amplification of 3 targeted sequences in the following genes :

- Alien gene : a non-natural sequence used as internal control for the amplification process – positions 1 and 4 on the strip of microtubes (see Figure 1).
- RdRP Gene : target specific to betacoronavirus (β CoV) human pathogens belonging to B and C groups (SARS-CoV-2, SARS-CoV and MERS-CoV) - positions 2 and 5 on the strip of microtubes (see Figure 1)
- N Gene : target specific to SARS-CoV-2 - positions 3 and 6 on the strip of microtubes (see Figure 1).

These amplifications are performed in parallel and the amplification result is obtained after a period of between 5 and 20 minutes depending on the viral load of the sample being analysed.

The BEC Reader device indicates in real time for each tube whether the amplification is positive (detection of the RNA target) or negative (no detection of the RNA target).

► **Operating Protocol**

▷ **Preparation of the samples**

The test is compatible with air and surface samples.

- Surface sampling guidelines:
Remove the swab from the package. When applying pressure with the swab onto the surface, move in at least two different directions while rotating the swab stick.
For surface, the recommended swab surface area is 25cm².
- Air sample guidelines:
Collect air sample with a Coriolis® Air sampler (Bertin Technologies).
Remove the swab from the package. Rub the inner surface of Coriolis® cone with the swab by applying a pressure while moving in at least two different directions and rotating the swab stick.

▷ **Reagent preparation**

- Open the aluminum pouch and check that the buffers or lyophilized pellets contained in the various tubes are well located at the bottom of the tubes and that no tubes are open. If this is the case, do not use the pouch.
If necessary, tap the tubes on the benchtop to make sure that drops of buffer or lyophilized pellets

do not remain in the cap.

- Leave the reagents at room temperature for 5 minutes before use. If the solution contained in the green tube is crystallized, warm the tube in your hands for a few seconds.
- Place the different tubes on the cardboard rack provided with the kit, respecting the color codes, see Figure 2. **Do not open the tubes at this stage.**


Figure 2: Rack that is provided for easy handling of the tubes.

▷ **BEC Reader device preparation**

- Switch on the BEC Reader device by pushing  button in front of the device during 3 seconds. The screen turns on and the reader performs a self-test. The cover opens.


Figure 3: BEC Reader device with location of ON/OFF switch button.


- Press on the screen and then press on TEST symbol


- Select the test "BEC_Environmental Assay Vx.x"


- Complete the RunID : this is the identifier of the run result file.
- Press "Auto 1-8" to automatically assign a name to the tubes. Go to the next screen by pressing ✓

- From this stage, the reader goes into preheating until reaching the set temperature of 65 °C.

- When the setpoint is reached, the reader informs the user.

▷ Test execution

CAUTION : perform the steps 1 through 5 with the 1st sample, and then repeat these steps with the 2nd sample. If only one sample is available, perform steps 1 through 6 and leave the wells 4 through 6 empty.

1) Step 1 : Lysis and RNA extraction

- a) Open the green tube (BEC SARS-CoV-2 Lysing Buffer (green)) in position 1 on the rack (see Figure 2).
- b) Immerse the swab in the solution and turn and stir the swab in the solution. Squeeze the swab against the inner side of the tube and discard the swab.
- c) Close the green tube, shake it vigorously during 10 seconds and place it in position 2 on the rack.

DO NOT leave the sample to stand before proceeding to Step 2 to avoid degradation of the mix that would interfere with the test.

2) Step 2 : Neutralization

- a) Open the blue tube (BEC SARS-CoV-2 Neutralization Buffer (blue)) and the green tube in position 2 on the rack.
- b) Withdraw the entirety of the liquid from the blue tube with the syringe fitted with the needle and dispense it in the green tube.
DO NOT pour directly without using the needle.
- c) Lay the syringe aside.
- d) Close the blue tube and the green tube. Shake the green tube vigorously during 10 seconds and place it on the rack in position 3.

3) Step 3 : Filtration

- a) Tap the green tube on the benchtop to lower the droplets that are in the cap. Open the green tube and the violet tube (BEC SARS-CoV-2 Transfer Tubes (violet)) in position 3.
- b) Draw as much liquid as possible from the green tube using the syringe fitted with the needle. Be careful to draw liquid from the green tube and not foam.
- c) Take out the needle from the syringe and replace it with the filter. Discard the needle in an appropriate waste bin.
- d) Insert the filter into the violet tube. Gently push the syringe plunger until it stops. By gently pushing the plunger, the flow of the filtrate is facilitated.
- e) Close the green tube and the violet tube. Place the violet tube in position 4 on the rack. Discard the green tube in an appropriate waste bin.

4) Step 4 : Transfer

- a) Open the violet tube and the orange tube (BEC SARS-CoV-2 Reaction Buffer (orange)) in position 4 on the rack.
- b) Aspirate a 100 μL volume from the violet tube with the violet fixed-volume pipette fitted with a filter-tip and transfer it to the orange tube.
- c) Close the violet tube and the orange tube. Shake vigorously the orange tube during 10 seconds and place it in position 5 on the rack.

CAUTION : Orange tube should be used in step 5 and the incubation (step 6) should be started within 30 minutes after preparation of orange tube.

5) Step 5 : Dispensing in the microtubes

- a) Open the orange tube in position 5.
 - b) Open the microtubes of the strip.
 - c) Use the orange pipette fitted with a filter-tip to aspirate a 25 μ L volume from the orange tube and dispense this volume in each microtube corresponding to the tested sample (see Figure 1), which means :
 - In tubes 1, 2 and 3 for sample 1,
 - In tubes 4, 5 and 6 for sample 2.
- Use a new filter tip for each microtube**
- d) Discard the filter-tip, close the orange tube and the microtubes.
 - e) Discard the green, blue, violet and orange tubes.

CAUTION : Repeat steps 1 through 5 with the second sample before proceeding to the next step (incubation).

When the strip of microtubes is filled, proceed immediately to step 6.

6) Step 6 : Incubation

CAUTION : Make sure the BEC reader has reached the temperature set point (65°C).

- a) Shake the microtubes until the lyophilized pellets are completely dissolved. Then tap the tubes on the benchtop to lower all the droplets to the bottom of the tube.
- b) Insert the microtube strip into the BEC Reader device, placing the black end tab on the left side of the instrument (see Figure 4). Leave the wells 7 and 8 empty.



Figure 4: How to insert microtube strip in the BEC Reader device. Red arrows showing position of the black end tab. Blue arrows showing empty well positions 7 and 8.

- c) Immediately close the lid. The BEC Reader device locks and starts the program automatically.
- d) A progress bar allows you to follow the progress of the test. The result for each well is displayed in real time on the screen. At the end of the amplification, the message "test completed" appears and all the results are displayed.

- e) The lid opens automatically at the end of the amplification.
- f) Discard immediately the microtubes **without opening the microtubes.**

► Interpretation of results

► Pictograms in positions 1 through 3 and 4 through 6

For each of the two samples (see Figure 1 for the position), the analysis of the results shall be performed pursuant to the following table :

Positive Control	RdRP gene	N gene	Interpretation of results
			Detection of viral RNA specific to the betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) and to the SARS-CoV-2 virus Sample is SARS-CoV-2 positive
			Detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) No detection of viral RNA specific to the SARS-CoV-2 Sample is presumptive positive for SARS-CoV-2 Repeat the test with a new sample.
			Detection of viral RNA specific to SARS-CoV-2 virus No detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) Sample is SARS-CoV-2 positive

Positive Control	RdRP gene	N gene	Interpretation of results
			No detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) and to SARS-CoV-2 virus Sample is SARS-CoV-2 negative
			Detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) and to SARS-CoV-2 virus No detection of positive control gene Sample is SARS-CoV-2 positive
			Detection viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) No detection of viral RNA specific to SARS-CoV-2 No detection of positive control gene Sample is presumptive positive for SARS-CoV-2 Repeat the test with a new sample.
			Detection of viral RNA specific to SARS-CoV-2 virus No detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) No detection of positiv control gene Sample is SARS-CoV-2 positive
			RT-LAMP inhibition or reagent failure. Repeat the test with a new sample.

Note : The « unclear »  pictogram may appear on one or several of the positions 1 through 6. In this case, the test is not valid as the amplification did not proceed properly. The test must be repeated with a new sample. In case the « unclear » pictogram appears on two successive runs, contact Technical Support.

▷ **Pictograms in positions 7 and 8**

Positions 7 and 8 always show the « unclear »  pictogram.

▶ **Optional: exporting results**

The BEC reader device automatically saves the test results into its memory up to 50 tests.

The test results can be exported as files in .CSV and .pdf format. Both files can be opened on a PC.


Pdf files contain run identification information (test name, RunID, date and time) and the results that are displayed on the screen.

To access detailed results (including fluorescence data), .CSV files must be opened with the 'T8 ISO Desktop' software.


- Connect the provided USB stick in the plug in front of the device

- Press the "Results" icon on the Home Screen. The list of test results is displayed.

- Select the test results and press the Export Results icon  to export results to the USB key.

Be sure the export is completed before removing the USB stick.

► Technical Support

Before you contact Bertin Bioreagent Technical Support, please gather the following information :

- Product Name
- Batch Number
- Serial Number of the instrument
- Error messages (if any)

Technical Support contact details :

- +33 (0)139 306 036
- tech@bertin-bioreagent.com

► Bibliography

Aebischer, A., Wernike, K., Hoffmann, B., & Beer, M. (2014). Rapid genome detection of Schmallenberg virus and bovine viral diarrhoea virus by use of isothermal amplification methods and high-speed real-time reverse transcriptase PCR. *Journal of Clinical Microbiology*, 52(6), 1883–1892. <https://doi.org/10.1128/JCM.00167-14>

Chander, Y., Koelbl, J., Puckett, J., Moser, M. J., Klingele, A. J., Liles, M. R., Carrias, A., Mead, D. A., & Schoenfeld, T. W. (2014). A novel thermostable polymerase for RNA and DNA loop-mediated isothermal amplification (LAMP). *Frontiers in Microbiology*, 5(AUG). <https://doi.org/10.3389/fmicb.2014.00395>

Francois, P., Tangomo, M., Hibbs, J., Bonetti, E. J., Boehme, C. C., Notomi, T., Perkins, M. D., & Schrenzel, J. (2011). Robustness of a loop-mediated isothermal amplification reaction for diagnostic applications. *FEMS Immunology and Medical Microbiology*, 62(1), 41–48. <https://doi.org/10.1111/j.1574-695X.2011.00785.x>

Li, J., & Macdonald, J. (2015). Advances in isothermal amplification: Novel strategies inspired by biological processes. In *Biosensors and Bioelectronics* (Vol. 64, pp. 196–211). Elsevier Ltd. <https://doi.org/10.1016/j.bios.2014.08.069>

Nagamine, K., Hase, T., & Notomi, T. (2002). Accelerated reaction by loop-mediated isothermal amplification using loop primers. *Molecular and Cellular Probes*, 16(3), 223–229. <https://doi.org/10.1006/mcpr.2002.0415>

Nie, K., Qi, S. xiang, Zhang, Y., Luo, L., Xie, Y., Yang, M. jie, Zhang, Y., Li, J., Shen, H., Li, Q., & Ma, X. jun. (2012). Evaluation of a Direct Reverse Transcription Loop-Mediated Isothermal Amplification Method without RNA Extraction for the Detection of Human Enterovirus 71 Subgenotype C4 in Nasopharyngeal Swab Specimens. *PLoS ONE*, 7(12). <https://doi.org/10.1371/journal.pone.0052486>

Notomi, T., Okayama, H., Masubuchi, H., Yonekawa, T., Watanabe, K., Amino, N., & Hase, T. (2000). Loop-mediated isothermal amplification of DNA. In *Nucleic Acids Research* (Vol. 28, Issue 12).

Yan, L., Zhou, J., Zheng, Y., Gamson, A. S., Roembke, B. T., Nakayama, S., & Sintim, H. O. (2014). Isothermal amplified detection of DNA and RNA. In *Molecular BioSystems* (Vol. 10, Issue 5, pp. 970–1003). Royal Society of Chemistry. <https://doi.org/10.1039/c3mb70304e>

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Evaluation of a Direct Reverse Transcription Loop-Mediated Isothermal Amplification Method without RNA Extraction for the Detection of Human Enterovirus 71 Subgenotype C4 in Nasopharyngeal Swab Specimens. *PLoS ONE*, 7(12). <https://doi.org/10.1371/journal.pone.0052486>

Notomi, T., Okayama, H., Masubuchi, H., Yonekawa, T., Watanabe, K., Amino, N., & Hase, T. (2000). Loop-mediated isothermal amplification of DNA. In *Nucleic Acids Research* (Vol. 28, Issue 12).

Yan, L., Zhou, J., Zheng, Y., Gamson, A. S., Roembke, B. T., Nakayama, S., & Sintim, H. O. (2014). Isothermal amplified detection of DNA and RNA. In *Molecular BioSystems* (Vol. 10, Issue 5, pp. 970–1003). Royal Society of Chemistry. <https://doi.org/10.1039/c3mb70304e>

French industrial company specialized in particular in Life Sciences and Biodefense fields. Bertin Technologies is the Legal Manufacturer of the tests and markets these tests.

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Leading company in the field of veterinary infectious disease diagnostics. Enalees developed the RNA extraction technology used for the tests and the point-of-care test model.

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French biotech company specialized in infectious disease in-vitro diagnostics, particularly in respiratory diseases. C4Diagnostics participated in the development of the tests and markets these tests.

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