



IL-4 (pig)

For laboratory research only. Not for human or veterinary diagnostic use.

Bertin Bioreagent also markets pre-analytical products, EIA kits, antibodies & biochemicals for:

- ▶ **Cardiology / Hypertension**
- ▶ **Diabetes / Obesity**
- ▶ **Endocrinology / Metabolism**
- ▶ **Inflammation**
- ▶ **Pharmacology**
- ▶ **Psychopharmacology**
- ▶ **Nitric Oxide**
- ▶ **Oncology / Apoptosis**
- ▶ **Oxidative injury**
- ▶ **Cell signaling**
- ▶ **Drug metabolism**

Do not hesitate to contact our after-sales services for further information at tech@bertin-bioreagent.com

IL-4 (pig)
Enzyme Immunoassay kit
#A05420.96 wells

For research laboratory use only
Not for human diagnostic use

This assay has been developed
& validated by Bertin Bioreagent

Fabriqué en France
Made in France



#A11420
Version: 0118

Table of contents

| | |
|--|----|
| ▶ PRECAUTION FOR USE..... | 6 |
| ▶ BACKGROUND..... | 7 |
| ▶ PRINCIPLE OF THE ASSAY | 7 |
| ▶ MATERIALS AND EQUIPMENT REQUIRED..... | 10 |
| ▶ SAMPLE COLLECTION AND PREPARATION..... | 11 |
| ▶ REAGENT PREPARATION | 12 |
| ▶ ASSAY PROCEDURE | 15 |
| ▶ DATA ANALYSIS..... | 21 |
| ▶ ACCEPTABLE RANGE..... | 22 |
| ▶ TYPICAL RESULTS..... | 23 |
| ▶ CHARACTERISTICS | 24 |
| ▶ TROUBLESHOOTING | 25 |
| ▶ BIBLIOGRAPHY | 26 |

96 wells**Storage: +4 °C****Expiry date: stated on the package**

This kit contains:

| Designation | Colour of cap | Item # | Quantity per kit | Form |
|---|-------------------------|----------------|------------------|-------------|
| IL-4 (pig) precoated 96-well Strip Plate | Blister with zip | A08420.1 ea | 1 | - |
| Low Streptavidin Poly_HRP Tracer | Green | A04417.100 dtn | 1 | Liquid |
| IL-4 (pig) Biotin- labelled Antibody | Red | A03420.100 dtn | 1 | Liquid |
| IL-4 (pig) Standard | Blue with red septum | A06420.1 ea | 2 | Lyophilised |
| Poly_HRP EIA Buffer | Grey / Blue | A07410.1 ea | 1 | Lyophilised |
| Wash Buffer | Silver | A17000.1 ea | 1 | Liquid |
| Tween 20 | Transparent | A12000.1 ea | 1 | Liquid |
| HRP Substrate Solution | Black | A09034.100 dtn | 1 | Liquid |
| HRP Stop Solution | Yellow | A22410.100 dtn | 1 | Liquid |
| Technical Booklet | - | A11420.1 ea | 1 | - |
| Well cover Sheet | - | - | 1 | - |

Each kit contains sufficient reagents for 96 wells. This allows for the construction of one standard curve in duplicate and the assay of 37 samples in duplicate.

▶ **Precaution for use**

Users are recommended to carefully read all instructions for use before starting work.

Each time a new pipette tip is used, aspirate a sample or reagent and expel it back into the same vessel. Repeat this operation two or three times before distribution in order to equilibrate the pipette tip.

- For research laboratory use only
- Not for human diagnostic use
- Do not pipet liquids by mouth
- Do not use kit components beyond the expiration date
- Do not eat, drink or smoke in area in which kit reagents are handled

▶ **Avoid splashing**

HRP Stop Solution and HRP Substrate Solution are harmful solutions. To avoid any contact, wear eye, hand, face and clothing protection when handling these.

Wearing gloves, laboratory coat and glasses is recommended when assaying kit materials and samples.

▶ **Temperature**

Unless otherwise specified, all the experiments are done at room temperature (RT), that is around +20°C. Working at +25°C or more affects the assay and decreases its efficiency.

► Background

Discovered in 1992, Interleukin-4 (IL-4) is a monomeric glycosylated polypeptide with a molecular weight of 18-20 kDa¹. IL-4 structure contains four alpha-helices which are anti-parallel and two anti-parallel beta-sheets²⁻³.

IL-4 is also known as B-cell stimulatory factor 1 or B-cell growth factor 1 and as Lymphocyte stimulatory factor.

IL-4 is a key regulator in humoral and adaptive immunity¹. IL-4 plays an essential role by promoting Th2 cell differentiation while inhibiting Th1 cell differentiation⁴. IL-4 is an anti-apoptosis factor for the immune and non-immune cells. Due to its main function, IL-4 is an anti-inflammatory cytokine in autoimmune diseases⁵. IL-4 deregulation could be a factor in the development of some diseases like allergic diseases⁴ or in the treatment resistance of tumors⁶.

► Principle of the assay

This Enzyme Immunometric Assay (EIA/ELISA) is based on a sandwich technique. The wells of the plate supplied are coated with a polyclonal antibody specific to IL-4 (pig).

IL-4 (pig) introduced into the wells (standard or sample) will be bound by the polyclonal antibody coated on the plate and is then detected by a second polyclonal antibody tagged with biotin also specific for IL-4 (pig).

The two antibodies then form a sandwich by binding on different parts of the IL-4 (pig).


The sandwich is immobilised on the plate so reagents in excess may be washed away.

The immunological complex is revealed by the interaction between biotin and streptavidin labelled with HRP (Tracer).

The concentration of IL-4 (pig) is determined by measuring the enzymatic activity of immobilized Tracer using TMB. The Tracer acts on TMB to form a yellow compound after the reaction has been stopped.

The intensity of the colour, which is determined by spectrophotometry at 450 nm, is proportional to the amount of IL-4 (pig) present in the well during the immunological incubation.

The principle of the assay is summarised below:



► **Materials and equipment required**

In addition to standard laboratory equipment, the following material is required:

For the assay:

- Precision micropipettes (20 to 1000 μL)
- Spectrophotometer plate reader (450 nm filter)
- Microplate washer (or washbottles)
- Orbital microplate shaker
- Multichannel pipette and disposable tips 30-300 μL
- UltraPure water #A07001.1L
- Polypropylene tubes



Water used to prepare all ELISA reagents and buffers must be UltraPure (deionized & free from organic contaminant traces).

Do not use distilled water, HPLC-grade water or sterile water.

- UltraPure water may be purchased from Bertin Bioreagent (item #A07001.1L).

▶ **Sample collection and preparation**

This assay may be used to measure IL-4 (pig).



It is the responsibility of the user to check the compatibility of the assay with the study matrix.

▶ **General precautions**

- All samples must be free from organic solvents prior to assay.
- Samples should be assayed immediately after collection or should be stored at -20°C or at -80°C .

▶ Reagent preparation

Each kit contains sufficient reagents for 96 wells. This allows for the construction of one standard curve in duplicate and the assay of 37 samples in duplicate.

All reagents need to be brought to room temperature (around +20°C) prior to the assay.

▶ Poly_HRP EIA Buffer

Reconstitute the Poly_HRP EIA Buffer #A07410 with 25 mL of UltraPure water. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion.

Stability at 4°C: 1 month.



Before use, filter the Buffer on 0,22 µm filter.

▶ IL-4 (pig) Standard

Reconstitute the IL-4 (pig) Standard vial #A06420 with 1 mL of UltraPure water. Allow it to stand for 5 minutes until completely dissolved and then mix thoroughly by gentle inversion.

The concentration of the first standard (S1) is 2.00 ng/mL. Prepare seven polypropylene tubes (for the seven other standards) and add 500 µL of Poly_HRP EIA Buffer into each tube. Then prepare the standards by serial dilutions as follow:

| Standard | Volume of Standard | Volume of Poly_HRP EIA Buffer | Standard concentration |
|----------|--------------------|-------------------------------|------------------------|
| S1 | - | - | 2.00 ng/mL |
| S2 | 500 μ L of S1 | 500 μ L | 1.00 ng/mL |
| S3 | 500 μ L of S2 | 500 μ L | 0.50 ng/mL |
| S4 | 500 μ L of S3 | 500 μ L | 0.25 ng/mL |
| S5 | 500 μ L of S4 | 500 μ L | 0.13 ng/mL |
| S6 | 500 μ L of S5 | 500 μ L | 0.06 ng/mL |
| S7 | 500 μ L of S6 | 500 μ L | 0.03 ng/mL |
| S8 | 500 μ L of S7 | 500 μ L | 0.02 ng/mL |

Stability at 4°C: within the day.

▶ **IL-4 (pig) Biotin-labelled Antibody**

The supplied IL-4 (pig) Biotin-labelled Antibody is concentrated 10 times. Calculate the volume needed (number of wells multiplied by 0.1 mL) and then dilute the IL-4 (pig) Biotin-labelled Antibody solution #A03420 with the appropriate volume of Poly_HRP EIA Buffer.

Example: for 40 wells you need 4 mL of IL-4 (pig) Biotin-labelled Antibody (40 x 0.1 mL), add 0.4 mL of IL-4 (pig) Biotin-labelled Antibody in 3.6 mL of Poly_HRP EIA Buffer.

Stability of diluted antibody at +4°C: within the day.

▶ **Wash Buffer**

Dilute 2 mL of concentrated Wash Buffer #A17000 with 800 mL of UltraPure water. Add 400 μ L of Tween 20 #A12000. Use a magnetic stirring bar to mix the content.

Stability at +4°C: 1 month.

▶ Assay procedure

It is recommended to perform the assays in duplicate following the instructions hereafter.

▶ Plate preparation

Prepare the Wash Buffer as indicated in the reagent preparation section.

Open the plate packet and select the sufficient strips for your assay and place the unused strips back in the packet.

Stability at +4°C: 1 month.

Rinse each well 5 times with Wash Buffer (300 μ L/well).

Just before distributing the reagents and samples, remove the buffer from the wells by inverting the plate and shaking out the last drops on a paper towel.

▶ Plate set-up

A plate set-up is suggested hereafter. The contents of each well may be recorded on the template sheet provided at the end of this technical booklet.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|----|----|---|---|---|---|---|---|----|----|----|
| A | Bk | S7 | S3 | * | * | * | * | * | * | * | * | * |
| B | Bk | S7 | S3 | * | * | * | * | * | * | * | * | * |
| C | Bk | S6 | S2 | * | * | * | * | * | * | * | * | * |
| D | NSB | S6 | S2 | * | * | * | * | * | * | * | * | * |
| E | NSB | S5 | S1 | * | * | * | * | * | * | * | * | * |
| F | NSB | S5 | S1 | * | * | * | * | * | * | * | * | * |
| G | S8 | S4 | * | * | * | * | * | * | * | * | * | * |
| H | S8 | S4 | * | * | * | * | * | * | * | * | * | * |

Bk : Blank

S1-S8 : Standards 1-8

NSB : Non Specific Binding

* : Samples

▶ Pipetting the reagents

All samples and reagents must reach room temperature prior to performing the assay.

Use different tips to pipet the buffers, standards, samples, antibody and other reagents.

Before pipetting, equilibrate the pipette tips in each reagent. Do not touch the liquid already in the well when expelling with the pipette tip.

> ***Poly_HRP EIA Buffer***

Dispense 100 μL to Non Specific Binding wells (NSB) wells.

> ***IL-4 (pig) Standard***

Dispense 100 μL of each of the eight standards (S8 to S1) in duplicate to appropriate wells.

Start with the lowest concentration standard (S8) and equilibrate the tip in the next higher standard before pipetting.

> ***IL-4 (pig) Sample***

Dispense 100 μL in duplicate to appropriate wells. Highly concentrated samples may be diluted in Poly_HRP EIA Buffer.

▷ **Incubating the plate**

Cover the plate with the cover sheet and incubate 60 minutes at room temperature, shaking at 300 rpm.

▷ **Washing the plate**

Rinse each well 5 times with Wash Buffer (300 μL /well). Just before distributing reagents, remove the buffer from the wells by inverting the plate and shaking out the last drops on a paper towel.

▷ **Pipetting the reagents**

> *IL-4 (pig) Biotin-labelled antibody*

Dispense 100 μ L to each well, except Blank (Bk) wells.

▷ **Incubating the plate**

Cover the plate with the cover sheet and incubate 60 minutes at room temperature, shaking at 300 rpm.

▷ **Washing the plate**

Rinse each well 5 times with Wash Buffer (300 μ L/well). Just before distributing reagents, remove the buffer from the wells by inverting the plate and shaking out the last drops on a paper towel.

▷ **Pipetting the reagents**

> *Low Streptavidin Poly_HRP Tracer*

Dispense 100 μ L to each well, except Blank (Bk) wells.

▷ **Incubating the plate**

Cover the plate with the cover sheet and incubate 30 minutes at room temperature, shaking at 300 rpm.

▷ **Developing and reading the plate**

- Empty the plate by turning it over. Rinse each well 5 times with 300 μ L of Wash Buffer. At the end of the last washing step, empty the plate and blot the plate on a paper towel to discard any trace of liquid.
- Add 100 μ L of HRP Substrate Solution to each well.
- Incubate the plate in the dark at room temperature without shaking. For the time, look at the lot specific Quality Control Sheet (QCS). In general, revelation time is 10 min.
- Add 100 μ L of HRP Stop Solution to each well.
- Wipe the bottom of the plate with a paper towel, and make sure that no liquid has splashed outside the wells.
- Read the plate at 450 nm (yellow color).

| Enzyme Immunoassay Protocol (volumes are in μL) | | | | |
|--|-------|-----|----------|--------|
| | Blank | NSB | Standard | Sample |
| Poly_HRP EIA Buffer | - | 100 | - | - |
| Standard | - | - | 100 | - |
| Sample | - | - | - | 100 |
| Cover plate, incubate 60 minutes at room temperature under orbital shaking at 300 rpm | | | | |
| Wash strips 5 times with 300 μL /well Discard liquid from the wells & dry on absorbent paper | | | | |
| Biotin-labelled Antibody | - | 100 | | |
| Cover plate, incubate 60 minutes at room temperature under orbital shaking at 300 rpm | | | | |
| Wash strips 5 times with 300 μL /well Discard liquid from the wells & dry on absorbent paper | | | | |
| Low Streptavidin Poly_HRP Tracer | - | 100 | | |
| Cover plate, incubate 30 minutes at room temperature under orbital shaking at 300 rpm | | | | |
| Wash strips 5 times with 300 μL /well Discard liquid from the wells & dry on absorbent paper | | | | |
| HRP Substrate Solution | 100 | | | |
| Incubate the plate in the dark without agitation | | | | |
| HRP Stop Solution | 100 | | | |
| Read the plate at 450 nm | | | | |

► Data analysis

Make sure that your plate reader has subtracted the absorbance readings of the blank wells from the absorbance readings of the rest of the plate. If not, do it now.

- Calculate the average absorbance for each NSB, standard and sample.
- For each standard, plot the absorbance on y axis versus the concentration on x axis. Draw a best-fit line through the points.
- To determine the concentration of your samples, find the absorbance value of each sample on the y axis.
- Read the corresponding value on the x axis which is the concentration of your unknown sample.
- Samples with a concentration greater than 2.00 ng/mL should be re-assayed after dilution in Poly_HRP EIA Buffer.
- Most plate readers are supplied with a curve-fitting software capable of graphing these data (4-parameter logistic fit 4PL). If you have this type of software, we recommend using it. Refer to it for further information.

► **Acceptable range**

- NSB absorbance ≤ 0.100 A.U.
- Limit of detection ≤ 0.02 ng/mL


► Typical results

The following data are for demonstration purpose only. Your data may be different and still correct.

These data were obtained using all reagents as supplied in this kit under the following conditions: 10 minutes developing at room temperature, reading at 450 nm. A 4 parameter logistic fitting with a 1/Y ponderation was used to determine the concentrations.

| Standard | IL-4 (pig) ng/mL | Absorbance A.U. |
|----------|------------------|-----------------|
| S1 | 2.00 | 1.847 |
| S2 | 1.00 | 1.245 |
| S3 | 0.50 | 0.695 |
| S4 | 0.25 | 0.380 |
| S5 | 0.13 | 0.225 |
| S6 | 0.06 | 0.132 |
| S7 | 0.03 | 0.086 |
| S8 | 0.02 | 0.063 |
| NSB | 0.00 | 0.050 |

Typical **IL-4 (pig)** standard curve



► Characteristics

- > **Limit of detection** calculated as the concentration of IL-4 (pig) corresponding to the NSB average plus three standard deviations is ≤ 0.02 ng/mL.

- > **Cross-reactivity**

| Species | Cross-reactivity |
|----------------------------|------------------|
| Recombinant IL-4 (bovine) | Moderate |
| Recombinant IL-4 (chicken) | None |
| Recombinant IL-4 (dolphin) | Strong |
| Recombinant IL-4 (equine) | None |
| Recombinant IL-4 (feline) | Moderate |
| Recombinant IL-4 (human) | None |
| Recombinant IL-4 (mouse) | Weak |
| Recombinant IL-4 (rabbit) | Weak |

► Troubleshooting

> ***Absorbance values are too low:***

- one reagent has not been dispensed,
- incorrect preparation,
- assay performed before reagents reached room temperature,
- reading time not long enough.

> ***High signal and background in all wells:***

- inefficient washing,
- overdeveloping (incubation time should be reduced),
- high ambient temperature.

> ***High dispersion of duplicates:***

- poor pipetting technique
- irregular plate washing.

These are a few examples of troubleshooting that may occur. If you need further explanation, Bertin Bioreagent will be happy to assist you. Feel free to contact our Technical Support staff by phone +33 (0)139 306 036, fax +33 (0)139 306 299 or E-mail tech@bertin-bioreagent.com, and be sure to indicate the batch number of the kit (see outside the box).

Bertin Bioreagent proposes EIA Training kit #B05005 and EIA workshop upon request. For further information, please contact our Technical Support.

► Bibliography

1. Stephen T. Smiley, Michael J. Grusby.
Interleukin 4.
1998, Encyclopedia of Immunology 1451–1453
2. Carr C, Aykent S, Kimack NM, Levine AD.
Disulfide assignments in recombinant mouse and human interleukin 4.
1991, Biochemistry. 30 (6): 1515–23
3. Walter MR, Cook WJ, Zhao BG *et al*.
Crystal structure of recombinant human interleukin-4.
1992, J. Biol. Chem. 267 (28): 20371–6
4. J. Zamorano, M.D. Rivas, M. Pérez-G.
Interleukin-4: A multifunctional cytokine.
2003, Inmunología Vol. 22/Num 2/215-224
5. Feldmann M, Brennan FM, Maini RN.
Role of cytokines in rheumatoid arthritis.
1996, Annu.Rev Immunol; 14: 397-440
6. Zamorano J, Kelly AE, Austrian J *et al*.
Costimulation of resting B lymphocytes alters the IL-4-activated IRS2 signaling pathway in a STAT6 independent manner: implications for cell survival and proliferation.
2011 Cell Res. 11: 44-54

| | | | | | | | | | | | | | | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|----|----|----|---|---|---|---|---|---|---|---|---|---|---|---|
| 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | |
| A | B | C | D | E | F | G | H | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ |



Bertin Bioreagent, over the last decades, has been developing and marketing over 100 biomarker assays, pre-analytical products, kits, antibodies and biochemicals thanks to its innovative work in research and development.

Our core areas are orientated to inflammation, oxidative injury, endocrinology, diabetes, obesity, hypertension, neurodegenerative diseases, HIV, prion diseases, pharmacokinetics and metabolism.

Bertin Bioreagent is active worldwide either with direct sales or through our qualified and trained international distribution network from the United States to Japan.

We are able to provide you with local technical support to use at ease our products.

For further information, please send your request to: tech@bertin-bioreagent.com



CONTACT US

Bertin Technologies
10 bis Avenue Ampère
Parc d'Activités du Pas du Lac
78180 Montigny-le-Bretonneux
FRANCE



+33 (0)139 306 036



tech@bertin-bioreagent.com



EU webstore: Bertin-bioreagent.com
US webstore: Bertin-corp.com